


ADAM-MC2, ADAM-CellIT

Automated fluorescence cell counter

ADAM-MC2, versatile automated fluorescence cell counter, allows users to perform assays for cells including cell counting and viability. To count cells, the cells are mixed with a Propidium Iodide (PI) stain and directly pipetted on to a disposable plastic chip. **ADAM-CellIT** with the same hardware as ADAM-MC2 complies 21 CFR part 11 which is a regulation about electronic records and signatures for use in cGMP facilities.


Principle of Viability Measurement (PI-Staining Method)



* There are two types of disposable chips: 2 channel and 4 channel

- After the samples are stained with fluorescent dye, propidium iodide, which intercalates DNA to stain the nucleus of target cells, ADAM-MC2 takes fluorescent images automatically.

QC Platform for producing CAR-T cell using ADAM-MC2, ADAM-CellIT



- For R&D, Process control, Quality control of CAR-T cell using ADAM-MC2, ADAM-CellIT

Specification

Specification	ADAM-MC2 & CellIT
Loading Vol	13 μ L (4ch)
Measurement Vol	3.4 μ L (4ch)
Optics	4 X
Analysis time	25 sec / test (4ch)
GMP & 21 CFR Part 11	ADAM CellIT (only)

ADAM CellIT ADAM MC2



ADAM-LS

Fluorescence cell analyzer with open tool system

ADAM-LS, versatile fluorescence cell analyzer based on 4-channel (bright field, UV LED, blue LED, green LED) for life science laboratory, allows users to perform lots of assays for cells including cell counting (total & nucleus), viability, apoptosis, cell cycle, fluorescence expression with open tool system and similarity to FACS.

Advantages

Versatility

- Compatible with a wide variety of eukaryotic cells
- Performs assays including cell counting (total & nucleus), viability, apoptosis, cell cycle, and fluorescence expression

High-quality cell analysis

- Performs 4-channel (bright field, UV, blue, and green LED)
- Result with histograms and select cell size gating

Open tool system

- Possibility to apply with other reagents in the right wavelength
- Fluorescence expression (GFP, RFP, DAPI)
- Nucleus counting (GFP, RFP, DAPI)

Similarity to FACS

- Dot plot
- Easy GUI, report format, and compatibility with FCS file
- No need for system maintenance



ADAM-LS specification

Optic channel	BR (Bright field) channel	FL (Fluorescence) channel 1	FL channel 2	FL channel 3
Channel	BR	DAPI	GFP	RFP
Light source	Green LED (up-light)	UV LED	Blue LED	Green LED
Excitation	-	FF01-390/40	FF01-466/40	FF01-543/22-25
Dichroic	FF495-Di03	Di02-R405	FF495-Di03	FF580-FDi01
Emission	FF01-525/50	FF01-452/45	FF01-525/50	BLP01-561R

ADAM-CDx

All-in-one system for cell therapy
(R&D and QC Process)

ADAM-CDx, all-in-one system with 4-channel (Bright field and 3 fluorescent channels) for cell therapy R&D and manufacturing process.


Advantages

Target cell

- Immune cell, CAR-T cell, NK cell
- Mammalian cell, Primary cell, Cell line
- Stem cells, PBMCs

High-quality All-in-one system

- Absolute counting of CD (cluster of differentiation) marker positive cells
- Total cell counting (BF, FL), Viability, Cell growth curve
- Apoptosis, Transfection efficiency, Cell cycling



ADAM-CDx specification

FL channel	BR (Bright Field) channel	FL (Fluorescence) channel 1	FL (Fluorescence) channel 2	FL (Fluorescence) channel 3
Dye	BR	BV510	PE	PerCP
Light source	Green LED (up-light)	405 (UV LED)	525 (Green)	525 (Green)
Excitation	-	390/40	510/42	525/50
Emission	572/28	525/50	572/28	BLP01-635R
CD marker panel (candidate)				
Panel 1 (T cell)		CD4	CD8	CD3
Panel 2 (T & B cell)		CD3	CD19	CD45
Panel 3 (NK cell)		CD3	CD56&16	CD45